(D)	Blutdruckmessgerät	
	Gebrauchsanweisung	2–19
(GB)	Blood pressure monitor	
	Instructions for use	20–35
(F)	Tensiomètre	
	Mode d'emploi	36–53
(E)	Tensiómetro	
	Instrucciones de uso	54–70
(I)	Misuratore di pressione	
	Istruzioni per l'uso	71–87
(TR)	Tansiyon ölçme cihazı	
	Kullanım kılavuzu	88–102
(RUS)	Прибор для измерения артериального давления на запястье	
	Инструкция по применению	103–121
(PL)	Ciśnieniomierz	
	Instrukcja obsługi	122–138
	Electromagnetic Compatibility Information	139–141

Inhalt

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Beurer GmbH
Servicecenter
Lessingstraße 10 b
89231 Neu-Ulm
Germany

Contents

1. Getting to know your device	20
2. Important notes	21
3. Unit description	24
4. Preparing for the measurement	25
5. Measuring blood pressure	27
6. Display and delete measured values	30
7. Error messages/troubleshooting	33
8. Cleaning and storing the device and cuff	34
9. Technical specifications	34

Included in delivery

- Blood pressure monitor with cuff
- 2 x 1.5 V LR03 AAA batteries
- Storage box
- Instructions for use

Dear customer,

Thank you for choosing one of our products. Our name stands for high-quality, thoroughly tested products for applications in the areas of heat, gentle therapy, blood pressure/diagnosis, weight, massage, beauty and air.

Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the notes they contain.

With kind regards, your Beurer team

1. Getting to know your device

Check that the packaging of the Beurer BC 57 blood pressure monitor has not been tampered with and make sure that all the required contents are present.

The wrist blood pressure monitor is used to carry out non-invasive measurement and monitoring of arterial blood pressure values in human adults. This allows you to quickly and easily measure your blood pressure, save the measurements and display the development of the measurements. You are also warned of possible existing cardiac arrhythmia.

The recorded values are classified and evaluated graphically.

2. Important notes

Signs and symbols

The following symbols are used in these instructions for use, on the packaging and on the type plate for the device and the accessories:

	Attention
	Note Note on important information
	Observe the instructions for use
	Application part, type BF
	Direct current
	Disposal in accordance with the Waste Electrical and Electronic Equipment EC Directive – WEEE
	Manufacturer

Storage 	Permissible storage temperature and humidity
Operating 	Permissible operating temperature and humidity
	Protect from moisture
	Serial number
	The CE labelling certifies that the product complies with the essential requirements of Directive 93/42/EEC on medical products.
	Certification symbol for products that are exported to the Russian Federation and members of the CIS

Notes on use

- In order to ensure comparable values, always measure your blood pressure at the same time of day.
- Before every measurement, relax for about five minutes.

- If you want to perform several measurements on the same person, wait five minutes between each measurement.
- Do not take a measurement within 30 minutes of eating, drinking, smoking or exercising.
- Repeat the measurement if you are unsure of the measured value.
- The measurements taken by you are for your information only – they are no substitute for a medical examination.

Discuss the measurements with your doctor and never base any medical decisions on them (e.g. medicines and their administration).

- Do not use the blood pressure monitor on newborns or patients with pre-eclampsia. We recommend consulting a doctor before using the blood pressure monitor during pregnancy.
- In the case of restricted circulation on the arm as a result of chronic or acute vascular diseases (including vascular constriction), the accuracy of the wrist measurement is limited. In this case you should avoid using an upper arm blood pressure monitor.
- Cardiovascular diseases may lead to incorrect measurements or have a detrimental effect on measurement accuracy. The same also applies to very low blood pressure, diabetes, circulatory disorders and arrhythmias as well as chills or shaking.

- This device is not intended for use by people (including children) with restricted physical, sensory or mental skills or a lack of experience and/or a lack of knowledge, unless they are supervised by a person who is responsible for their safety or are instructed by such a person in how to use the device. Supervise children around the device to ensure they do not play with it.
- The blood pressure monitor must not be used in connection with a high-frequency surgical unit.
- Only use the unit on people who have the specified wrist measurement for the device.
- Please note that when inflating, the functions of the limb in question may be impaired.
- During the blood pressure measurement, the blood circulation must not be stopped for an unnecessarily long time. If the device malfunctions, remove the cuff from the arm.
- Do not allow sustained pressure in the cuff or frequent measurements. The resulting restriction of the blood flow may cause injury.
- Ensure that the cuff is not placed on an arm in which the arteries or veins are undergoing medical treatment, e.g. intravascular access or therapy, or an arterio-venous (AV) shunt.
- Do not use the cuff on people who have undergone a mastectomy.

- Do not place the cuff over wounds as this may cause further injury.
- Place the cuff on your wrist only. Do not place the cuff on other parts of the body.
- The blood pressure monitor can only be operated with batteries.
- To conserve the batteries, the blood pressure monitor switches off automatically if you do not press any buttons for 30 seconds.
- The device is only intended for the purpose described in these instructions for use. The manufacturer is not liable for damage resulting from improper or careless use.
- Please note that data transfer and data storage is only possible when your blood pressure monitor is supplied with power. As soon as the batteries are empty, the blood pressure monitor loses the date and time.

Instructions for storage and maintenance

- The blood pressure monitor is made from precision and electronic components. The accuracy of the measurements and service life of the device depend on its careful handling:
 - Protect the device from impacts, humidity, dirt, marked temperature fluctuations and direct sunlight.
 - Do not drop the device.
- Do not use the device in the vicinity of strong electromagnetic fields and keep it away from radio systems or mobile telephones.
- We recommend that the batteries be removed if the device will not be used for a prolonged period of time.

Notes on handling batteries

- If your skin or eyes come into contact with battery fluid, flush out the affected areas with water and seek medical assistance.
- **Choking hazard!** Small children may swallow and choke on batteries. Store the batteries out of the reach of small children.
- Observe the plus (+) and minus (-) polarity signs.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.
- Protect the batteries from excessive heat.
- **Risk of explosion!** Never throw batteries into a fire.
- Do not charge or short-circuit batteries.
- If the device is not to be used for a long period, take the batteries out of the battery compartment.
- Use identical or equivalent battery types only.
- Always replace all batteries at the same time.
- Do not use rechargeable batteries.
- Do not disassemble, split or crush the batteries.


Instructions for repairs and disposal

- Batteries do not belong in household waste. Please dispose of empty batteries at the collection points intended for this purpose.
- Do not open the device. Failure to comply will result in voiding of the warranty.
- Do not repair or adjust the device yourself. Proper operation can no longer be guaranteed in this case.
- Repairs must only be carried out by Customer Services or authorised suppliers. Before making a claim, please check the batteries first and replace them if necessary.
- For environmental reasons, do not dispose of the device in the household waste at the end of its useful life. Dispose of the unit at a suitable local collection or recycling point. Dispose of the device in accordance with EC Directive – WEEE (Waste Electrical and Electronic Equipment). If you have any questions, please contact the local authorities responsible for waste disposal.


3. Unit description

1. Display
2. **START/STOP** button ①
3. Memory button **M**
4. Wrist cuff
5. Battery compartment lid
6. Risk indicator

Information on the display:

1. Risk indicator
2. Symbol for *Bluetooth*® transfer
3. Time and date
4. Systolic pressure
5. Diastolic pressure
6. Calculated pulse value
7. Cardiac arrhythmia symbol
- Pulse symbol
8. Battery status indicator
9. Release air (arrow)
10. User memory
11. Memory space number/memory display
average value (A), morning (AM), evening (PM)

System requirements for the PC software „beurer HealthManager“ :

- from Windows 7 SP1
- from USB 2.0 (Type A)

System requirements for the app “beurer Health-Manager”:

- iOS ≥ 8.0, Android™ ≥ 4.4
- *Bluetooth*® ≥ 4.0


List of compatible devices:

4. Preparing for the measurement

Inserting the batteries

- Remove the battery compartment lid on the left side of the device.
- Insert two 1.5 V micro (alkaline type LR03) batteries. Make sure that the batteries are inserted the correct way round. Do not use rechargeable batteries.
- Close the battery compartment lid again carefully.

If the battery change symbol  is flashing and **Lo** appears, no further measurements are possible and you must replace all batteries.

24h flashes on the display. Now set the date and time as described below.

The empty, completely flat batteries must be disposed of through specially designated collection boxes, recycling points or electronics retailers. You are legally required to dispose of the batteries.

The codes below are printed on batteries containing harmful substances:

Pb = Battery contains lead,

Cd = Battery contains cadmium,

Hg = Battery contains mercury.

Set the hour format, date, time and Bluetooth® settings

This menu allows you to set the following functions, one after another.

It is essential to set the date and time. Otherwise, you will not be able to save your measured values correctly with a date and time and access them again later.

i If you press and hold the memory button **M**, you can set the values more quickly.

Press and hold the **START/STOP** **①** button for 5 seconds if you have not reinserted the batteries.

Hour format

The hour format now flashes on the display.

- Select the desired hour format using the memory button **M** and confirm with the **START/STOP** button **①**.

Date

The year flashes on the display.

- Select the desired year using the memory button **M** and confirm with the **START/STOP** button **①**.

The month flashes on the display.

- Select the desired month using the memory button **M** and confirm with the **START/STOP** button **①**.

The day flashes on the display.

- Select the desired day using the memory button **M** and confirm with the **START/STOP** button **①**.

i If the hour format is set as **12h**, the day/month display sequence is reversed.

The hours flash on the display.

- Select the desired hour using the memory button **M** and confirm with the **START/STOP** button ①.

The minutes flash on the display.

- Select the desired minute using the memory button **M** and confirm with the **START/STOP** button ①.


The *Bluetooth*® symbol flashes on the display.

- Use the memory button **M** to select whether automatic *Bluetooth*® data transfer is activated (*Bluetooth*® symbol flashes) or deactivated (*Bluetooth*® symbol is not shown) and confirm with the **START/STOP** button ①.

i *Bluetooth*® transfers will reduce the battery life.

5. Measuring blood pressure

Attaching the cuff

- Expose your left wrist. Ensure that the circulation of the arm is not hindered by tight clothing or similar. Place the cuff on the inside of your wrist.
- Close the cuff with the hook-and-loop fastener such that the upper edge of the device is positioned approx. 1 cm underneath the heel of the hand.
- The cuff must be fitted tightly around the wrist, but must not constrict.

Adopt the correct posture

- Before every measurement, relax for about five minutes. Otherwise deviations can occur.
- You can take the measurement while sitting or lying. To take your blood pressure, make sure you are sitting comfortably with your arms and back leaning on something. Do not cross your legs. Place your

feet flat on the ground. Make sure to rest your arm and move it. Always make sure that the cuff is at heart level. Otherwise significant deviations can occur. Relax your arm and the palm of your hand.

- To avoid falsifying the measurement, it is important to remain still during the measurement and not to speak.

Performing the blood pressure measurement

As described above, attach the cuff and adopt the posture in which you want to perform the measurement.

- Press the **START/STOP** button ① to start the blood pressure monitor. All displays will illuminate briefly.

The blood pressure monitor will begin the measurement automatically after 5 seconds.


The measurement is taken during the inflation phase.

- ① You can cancel the measurement at any time by pressing the **START/STOP** button ①.

As soon as a pulse is found, the pulse symbol ♥ is displayed.

- The systolic pressure, diastolic pressure and pulse rate measurements are displayed.

- **EE_** appears if the measurement could not be performed properly. Observe the chapter on error messages/troubleshooting in these instructions for use and repeat the measurement.
- Now select the desired user memory by pressing the memory button **M**. If you do not select a user memory, the measurement is stored in the most recently used user memory. The relevant symbol **M1** or **M2** appears in the display.
- Press the **START/STOP** button ① to switch off the blood pressure monitor. The measurement is then stored in the selected user memory.

If **Bluetooth® data transfer** is activated, data is transferred after pressing the **START/STOP** button ①.

- The **Bluetooth®** symbol flashes on the display. The blood pressure monitor now attempts to connect to the app for approx. 30 seconds.

- The *Bluetooth*® symbol stops flashing as soon as a connection is established. All measurement data is automatically transferred to the app. Once the data has been successfully transferred, the device switches off.
- If a connection to the app cannot be established after 30 seconds, the *Bluetooth*® symbol goes out and the blood pressure monitor switches off automatically after 3 minutes.
- If you forget to turn off the device, it will switch off automatically after a maximum of 3 minutes. In this case too, the value is stored in the selected or most recently used user memory.
- Wait at least 5 minutes before taking another measurement.

Evaluating results

Cardiac arrhythmia:

This device can identify potential disruptions of the heart rhythm when measuring and if necessary, indicates this after the measurement with the symbol .

This can be an indicator for arrhythmia. Arrhythmia is an illness in which the heart rhythm is abnormal because of flaws in the bioelectrical system that regulates the heart-beat. The symptoms (skipped or premature heart beats, pulse being slow or too fast) can be caused by factors such as heart disease, age, physical make-up, excess

stimulants, stress or lack of sleep. Arrhythmia can only be determined through an examination by your doctor.

If the symbol  is shown on the display after the measurement has been taken, repeat the measurement. Please ensure that you rest for 5 minutes beforehand and do not speak or move during the measurement. If the symbol  appears frequently, please consult your doctor. Self-diagnosis and treatment based on the measurements can be dangerous. Always follow your GP's instructions.

Risk indicator:

The measurements can be classified and evaluated in accordance with the following table.

However, these standard values serve only as a general guideline, as the individual blood pressure varies in different people and different age groups etc.

It is important to consult your doctor regularly for advice. Your doctor will tell you your individual values for normal blood pressure as well as the value above which your blood pressure is classified as dangerous.

The bar chart on the display and the scale on the unit show which category the recorded blood pressure values fall into. If the values of systole and diastole fall into two different categories (e.g. systole in the 'High normal' category and diastole in the 'Normal' category), the graphical classification on the device always shows the higher category; for the example given this would be 'High normal'.

Blood pressure value category	Systole (in mmHg)	Diastole (in mmHg)	Action
Setting 3: severe hypertension	≥ 180	≥ 110	seek medical attention
Setting 2: moderate hypertension	160–179	100–109	seek medical attention
Setting 1: mild hypertension	140–159	90–99	regular monitoring by doctor
High normal	130–139	85–89	regular monitoring by doctor
Normal	120–129	80–84	self-monitoring
Optimal	< 120	< 80	self-monitoring

Source: WHO, 1999

6. Display and delete measured values

The results of every successful measurement are stored together with the date and time. The oldest measurement is overwritten in the event of more than 60 measurements.

- Once the device is switched off, select the desired user memory (1 2) using the memory button **M** and confirm with the **START/STOP** button ①. The average value of all measurements appears on the display.

If *Bluetooth®* is activated (the  symbol flashes on the display), the blood pressure monitor attempts to connect to the app. If you press the memory button **M** in the meantime, the transfer is cancelled and the average values are displayed. The  symbol is no longer shown. The buttons are deactivated and the  symbol ceases to flash as soon as a connection is established and the data is transferred.

User memory

A flashes on the display.

The average value of all saved measured values in this user memory is displayed.

- Press the memory button **M**.

AM flashes on the display.

The average value of the morning measurements for the last 7 days is displayed (morning: 5.00 a.m. – 9.00 a.m.).


- Press the memory button **M**.

PM flashes on the display.

The average value of the evening measurements for the last 7 days is displayed (evening: 6.00 p.m. – 8.00 p.m.).

- When the memory button **M** is pressed again, the last individual measurement is displayed (in this example, measurement 03).

- When the memory button **M** is pressed again, you can view your individual measured values.
- To switch the device off again, press the **START/STOP** button ①.

- ① You can exit the menu at any time by pressing the **START/STOP** button ①.

- To clear the memory of the relevant user memory, you must first select a user memory.
- Start the retrieval of the average measured values. **A** flashes on the display and the average value of all saved measured values in this user memory is displayed.
- Press and hold the memory button **M** for 5 seconds depending on the user memory you are in.

All the values in the current user memory are deleted.

To clear individual measurements from the relevant user memory, you must first select a user memory.

- Start individual measurement access.
- Press and hold the memory button **M** for 5 seconds (depending on the user memory you are in).
- The selected value is deleted. The device briefly displays **LL 00**.
- If you would like to delete other values, repeat the process described above.

You can switch off the device at any time by pressing the **START/STOP** button **1**.

Transfer via *Bluetooth*®

It is also possible to transfer the measured values saved on the device to your smartphone via *Bluetooth*®.

You will need the “beurer HealthManager” app for this. This is available in the Apple App Store and from Google Play.

Proceed as follows to transfer values:

If *Bluetooth*® is activated in the settings menu, the data is transferred automatically after the measurement. The **Bluetooth** symbol appears in the top left of the display (see chapter 4 “Preparing the measurement”).

Step 1: BC 57

Activate *Bluetooth*® on your device (see chapter 4 “Preparing the measurement, *Bluetooth*®”).

Step 2: “beurer HealthManager” app

In the “beurer HealthManager” app, add the BC 57 under “Settings/My devices”.

Step 3: BC 57

take a measurement.

Step 4: BC 57:

data transfer immediately following measurements.

- If **Bluetooth® data transfer** has been activated, data is transferred after having confirmed the user memory and pressing the **START/STOP** button ①.

Step 4: BC 57:

Data transfer at a later point:

- Go to memory mode (chapter 7). Select the desired user memory. The **Bluetooth®** transfer starts automatically.

- ❗ The “beurer HealthManager” app must be active to allow data transfer.

If your smartphone has a protective cover, remove this to ensure that there is no interference during the transfer. Begin the data transfer in the “beurer HealthManager” app.

We would explicitly draw attention to the fact that the software to hand is not a medical product in accordance with EU Directive 93/42/EEC.

The values displayed are purely for visualisation purposes and must not be used as a basis for therapeutic treat-

ments. The software is not part of a diagnostic-medical system.

7. Error messages/troubleshooting

In case of faults, the **EE_** message appears in the display. Error messages may appear if:

- systolic or diastolic pressure could not be measured (**EE 1** appears on the display)
- systolic or diastolic pressure was outside the measurement range (**H 1** or **L 0** appears on the display)
- the cuff is fastened too tightly or loosely (**EE 2** appears on the display)
- the pump pressure is higher than 300 mmHg (**EE 8** appears on the display)
- if the pressure during pumping up remains below 15 mmHg for longer than 3 minutes (**EE 3** appears on the display)
- there is a system or device error (**EE 0**, **EE 4**, **EE 5** or **EE 6** appears on the display)
- the data could not be sent via **Bluetooth®** (**EE 7**).

In such cases, repeat the measurement. Ensure that you do not move or speak. If necessary, reinsert or replace the batteries.

8. Cleaning and storing the device and cuff

- Clean the device and cuff carefully using a slightly damp cloth only.
- Do not use any cleaning agents or solvents.
- Under no circumstances hold the device and cuff under water, as this can cause liquid to enter and damage the device and cuff.
- If you store the device and cuff, do not place heavy objects on the device and cuff. Remove the batteries.

9. Technical specifications

Model no.	BC 57
Measurement method	Oscillometric, non-invasive blood pressure measurement on the wrist
Measurement range	Cuff pressure 0-300 mmHg, Systolic 60-260 mmHg, Diastolic 40-199 mmHg, Pulse 40-180 beats/minute
Display accuracy	Systolic ± 3 mmHg, Diastolic ± 3 mmHg, Pulse ± 5 % of the value shown
Measurement inaccuracy	Max. permissible standard deviation according to clinical testing: Systolic 8 mmHg / Diastolic 8 mmHg

Memory	2 x 60 memory spaces
Dimensions	L 78 mm x W 60 mm x H 23 mm
Weight	Approximately 96 g (without batteries)
Cuff size	140 to 195 mm
Permissible operating conditions	+10 °C to +40 °C, ≤ 85 % relative humidity (non-condensing)
Permissible storage conditions	-20 °C to +50 °C, ≤ 85 % relative humidity, 800-1050 hPa ambient pressure
Power supply	2 x 1.5V --- AAA batteries
Battery life	For approx. 200 measurements, depending on levels of blood pressure and pump pressure
Accessories	Instructions for use, 2 x 1.5V AAA batteries, storage box
Data transfer	2402 MHz – 2480 MHz frequency band Transmission power max. 2.8 dBm The blood pressure monitor uses <i>Bluetooth</i> ® low energy technology Compatible with <i>Bluetooth</i> ® 4.0 smartphones/tablets List of supported smartphones/tablets

The serial number is located on the device or in the battery compartment.

- This device is in line with European Standard EN 60601-1-2 and is subject to particular precautions with regard to electromagnetic compatibility (EMC). Please note that portable and mobile HF communication systems may interfere with this unit. More details can be requested from the stated Customer Services address or found at the end of the instructions for use.
- This device corresponds to the EU Medical Devices Directive 93/42/EEC, the German Medical Devices Act (Medizinproduktegesetz) and the standards EN 1060-1 (non-invasive sphygmomanometers, Part 1: General requirements), EN 1060-3 (non-invasive sphygmomanometers, Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems) and IEC 80601-2-30 (Medical electrical equipment – Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers).
- The accuracy of this blood pressure monitor has been carefully checked and developed with regard to a long useful life.

If using the device for commercial medical purposes, it must be regularly tested for accuracy by appropriate

means. Precise instructions for checking accuracy may be requested from the service address.

- We hereby guarantee that this product complies with the European RED Directive 2014/53/EU. Please contact kd@beurer.de to obtain detailed information, such as the CE Declaration of Conformity.

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Electromagnetic Compatibility Information

For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacture's declaration – electromagnetic emission

The BC 57 is intended for use in the electromagnetic environment specified below.
The customer or the user of the BC 57 should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 2	The BC 57 must emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected. The BC 57 is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacture's declaration – electromagnetic immunity

The BC 57 is intended for use in the electromagnetic environment specified below. The customer or the user of the BC 57 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

Guidance and manufacture's declaration – electromagnetic immunity

The BC 57 is intended for use in the electromagnetic environment specified below. The customer or the user of the BC 57 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the BC 57, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2.3 \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength

Guidance and manufacture's declaration – electromagnetic immunity

in the location in which the BC 57 is used exceeds the applicable RF compliance level above, the BC 57 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the BC 57.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and the BC 57

The BC 57 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the BC 57 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the BC 57 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

The *Bluetooth*® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners.

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries.

Google Play and the Google Play logo are trademarks of Google Inc.

Android is a trademark of Google Inc.

Beurer GmbH • Söflinger Straße 218 • 89077 Ulm, Germany • www.beurer.com
www.beurer-gesundheitsratgeber.com • www.beurer-healthguide.com