DEUTSCH

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ENGLISH

Signs and symbols

The following symbols are used in these instructions for use, on the packaging and on the type plate for the device and accessories:

	Caution		Follow instructions for use
	Permissible storage and transport temperature and humidity		Type BF applied part
	Permissible operating temperature and humidity		Disposal in accordance with EC Directive – WEEE (Waste Electrical and Electronic Equipment).
	Protected against ingress of solid foreign objects greater than 12.5 mm in diameter. Protected against drops of water falling at up to 15° from vertical.		The CE labelling certifies that the product complies with the essential requirements of Directive 93/42/EEC on medical products.
	The CE labelling certifies that the product complies with the essential requirements of Directive 93/42/EEC on medical products.		

Included in delivery

- Thermometer
- These instructions for use
- 1 x 1,5 V LR 41 battery

Important safety information – retain for later use.

In order to make optimum use of all the merits of the clinical thermometer, you should carefully read the operating instructions before using it, retain the instructions for later use and also keep them accessible for other users.

• The thermometer is intended to measure the temperature of the human body.

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→

ESPAÑOL

Guidance and manufacturer's declaration – electromagnetic emissions

The FT 15/1 is intended for use in the electromagnetic environment specified below. The customer or the user of the FT 15/1 should assure that it is used in such an environment.

Emissions test

RF emissions CISPR 11 Group 1

The FT 15/1 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

RF emissions CISPR 11 Class B

The FT 15/1 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Harmonic emissions IEC 61000-3-2 Not applicable

Voltage fluctuations/flicker emissions IEC 61000-3-3 Not applicable

→

Electromagnetic Compatibility Information

Table 1

Guidance and manufacturer's declaration – electromagnetic emissions

The FT 15/1 is intended for use in the electromagnetic environment specified below. The customer or the user of the FT 15/1 should assure that it is used in such an environment.

Immunity test

Compliance

Electromagnetic environment – guidance

The FT 15/1 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

The FT 15/1 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

If field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast can be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the FT 15/1 is used exceeds the applicable RF compliance level above, the FT 15/1 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the FT 15/1.

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